



You are currently viewing the documentation for the Opentrons OT-1. To view documentation for the OT-2, click [here](#).

API Reference

If you are reading this, you are probably looking for an in-depth explanation of API classes and methods to fully master your protocol development skills.

Robot

All protocols are set up, simulated and executed using a Robot class.

class `opentrons.Robot`

This class is the main interface to the robot.

Through this class you can can:

- define your `opentrons.Deck`
- `connect()` to Opentrons physical robot
- `home()` axis, move head (`move_to()`)
- `pause()` and `resume()` the protocol run
- set the `head_speed()` of the robot

Each Opentrons protocol is a Python script. When evaluated the script creates an execution plan which is stored as a list of commands in Robot's command queue.

Here are the typical steps of writing the protocol:

- Using a Python script and the Opentrons API load your containers and define instruments (see `Pipette`).

- Call `reset()` to reset the robot's state and clear commands.

- Write your instructions which will get converted into an execution plan.

- Review the list of commands generated by a protocol `commands()`.

- `connect()` to the robot and call `run()` it on a real robot.

See `Pipette` for the list of supported instructions.

Examples

```
>>> from opentrons import robot, instruments, containers
>>> robot.reset()
<opentrons.robot.Robot object at ...>
>>> plate = containers.load('96-flat', 'A1', 'plate')
>>> p200 = instruments.Pipette(axis='b', max_volume=200)
>>> p200.aspirate(200, plate[0])
<opentrons.instruments.Pipette object at ...>
>>> robot.commands()
['Aspirating 200 at <Deck><Slot A1><Container plate><Well A1>']
```

actions

Return a copy of a raw list of commands in the Robot's queue.



attached to.
 instrument (*Instrument*) – An instance of a **Pipette** to
 attached to the axis.

Notes

A canonical way to add to add a Pipette to a robot is:

```
from opentrons.instruments.pipette import Pipette
p200 = Pipette(axis='a')
```

This will create a pipette and call `add_instrument()` to attach the instrument.

commands()

Access the human-readable list of commands in the robot's queue.

Returns: *A list of string values for each command in the queue, for example*
 'Aspirating 200uL at <Deck>/<Slot A1>/<Container
 plate>/<Well A1>'

connect(port=None, options=None)

Connects the robot to a serial port.

Parameters: port (*str*) – OS-specific port name or 'Virtual
 Smoothie'
 options (*dict*) – if **port** is set to 'Virtual Smoothie',
 provide the list of options to be passed to
get_virtual_device()

Returns:
 Return **True** for success, **False** for failure.
 type:

containers()

Returns the dict with all of the containers on the deck.

diagnostics()

Access diagnostics information for the robot.

Returns: **axis_homed** — axis that are currently in home position.
switches — end stop switches currently hit.
steps_per_mm — steps per millimeter calibration
 values for **x** and **y** axis.

Return type: Dictionary with the following keys

disconnect()

Disconnects from the robot.

get_mosfet(mosfet_index)

Get MOSFET for a MagBead (URL).

Parameters: mosfet_index (*int*) – Number of a MOSFET on MagBead.
 Returns:
 Return type: Instance of **InstrumentMosfet**.

get_motor(axis)

**get_warnings()**

Get current runtime warnings.

Returns: Runtime warnings accumulated since the last `run()` or `simulate()`.

head_speed(*args, **kwargs)

Set the XY axis speeds of the robot, set in millimeters per minute

Parameters: rate (*int*) – An integer setting the mm/minute rate of the X and Y axis. Speeds too fast (around 6000 and higher) will cause the robot to skip step, be careful when using this method

Examples

```
>>> from opentrons import robot
>>> robot.connect('Virtual Smoothie')
>>> robot.home()
>>> robot.head_speed(4500)
>>> robot.move_head(x=200, y=200)
```

home(*args, **kwargs)

Home robot's head and plunger motors.

Parameters: *args –
A string with axes to home. For example 'xyz' or 'ab'.

If no arguments provided home Z-axis then X, Y, B, A

Notes

Sometimes while executing a long protocol, a robot might accumulate precision error and it is recommended to home it. In this scenario, add `robot.home('xyzab')` into your script.

Examples

```
>>> from opentrons import Robot
>>> robot.connect('Virtual Smoothie')
>>> robot.home()
```

move_to(location, instrument=None, strategy='arc', **kwargs)

Move an instrument to a coordinate, container or a coordinate within a container.

Parameters: location (*one of the following:*) – 1. **Placeable** (i.e. Container, Deck, Slot, Well) — will move to the origin of a container. 2. **Vector** move to the given coordinate in Deck coordinate system. 3. (**Placeable**, **Vector**) move to a given coordinate within object's coordinate system.
instrument – Instrument to move relative to. If `None`, move relative to the center of a gantry.
strategy (`{'arc', 'direct'}`) –
`arc` : move to the point using arc trajectory avoiding



Examples

```
>>> from opentrons import Robot
>>> robot.reset()
<opentrons.robot.robot.Robot object at ...>
>>> robot.connect('Virtual Smoothie')
>>> robot.home()
>>> plate = robot.add_container('96-flat', 'A1', 'plate')
>>> robot.move_to(plate[0])
>>> robot.move_to(plate[0].top())
```

pause()

Pauses execution of the protocol. Use [resume\(\)](#) to resume

reset()

Resets the state of the robot and clears:

- Deck
- Instruments
- Command queue
- Runtime warnings

resume()

Resume execution of the protocol after [pause\(\)](#)

stop()

Stops execution of the protocol.

Pipette

`class opentrons.instruments.Pipette(robot, axis, name=None, channels=1, min_volume=0, max_volume=None, trash_container=None, tip_racks=[], aspirate_speed=300, dispense_speed=500)`

Through this class you can:

- Handle liquids with [aspirate\(\)](#), [dispense\(\)](#), [mix\(\)](#), and [blow_out\(\)](#)
- Handle tips with [pick_up_tip\(\)](#), [drop_tip\(\)](#), and [return_tip\(\)](#)
- Calibrate this pipette's plunger positions
- Calibrate the position of each [Container](#) on deck

Here are the typical steps of using the Pipette:

Instantiate a pipette with a maximum volume (uL)

and an axis (*a* or *b*) * Design your protocol through the pipette's liquid-handling commands

Parameters:

- `axis (str)` – The axis of the pipette's actuator on the Opentrons robot ('a' or 'b')
- `name (str)` – Assigns the pipette a unique name for saving it's calibrations
- `channels (int)` – The number of channels on this pipette (Default: 1)
- `min_volume (int)` – The smallest recommended uL volume for this pipette (Default: 0)
- `max_volume (int)` – The largest uL volume for this pipette (Default: `min_volume + 1`)
- `trash_container (Container)` – Sets the default location [drop_tip\(\)](#) will put tips (Default: `None`)
- `tip_racks (list)` – A list of Containers for this Pipette to track tips when calling [pick_up_tip\(\)](#) (Default: [])
- `aspirate_speed (int)` – The speed (in mm/minute) the



Return type: A new instance of [Pipette](#).

Examples

```
>>> from opentrons import instruments, containers
>>> p1000 = instruments.Pipette(axis='a', max_volume=1000)
>>> tip_rack_200ul = containers.load('tiprack-200ul', 'A1')
>>> p200 = instruments.Pipette(
...     name='p200',
...     axis='b',
...     max_volume=200,
...     tip_racks=[tip_rack_200ul])
```

aspirate(*volume=None, location=None, rate=1.0*)

Aspirate a volume of liquid (in microliters/uL) using this pipette

Notes

If no *location* is passed, the pipette will aspirate from it's current position. If no *volume* is passed, *aspirate* will default to it's *max_volume*

Parameters: *volume* (*int or float*) – The number of microliters to aspirate (Default: self.max_volume)
location ([Placeable](#) or tuple([Placeable](#), [Vector](#))) – The [Placeable](#) ([Well](#)) to perform the aspirate. Can also be a tuple with first item [Placeable](#), second item relative [Vector](#)
rate (*float*) – Set plunger speed for this aspirate, where speed = rate * aspirate_speed (see [set_speed\(\)](#))

Returns:

Return type: This instance of [Pipette](#).

Examples

```
>>> p200 = instruments.Pipette(
...     name='p200', axis='a', max_volume=200)

>>> # aspirate 50uL from a Well
>>> p200.aspirate(50, plate[0])
<opentrons.instruments.pipette.Pipette object at ...>

>>> # aspirate 50uL from the center of a well
>>> p200.aspirate(50, plate[1].bottom())
<opentrons.instruments.pipette.Pipette object at ...>

>>> # aspirate 20uL in place, twice as fast
>>> p200.aspirate(20, rate=2.0)
<opentrons.instruments.pipette.Pipette object at ...>

>>> # aspirate the pipette's remaining volume (80uL) from a Well
>>> p200.aspirate(plate[2])
<opentrons.instruments.pipette.Pipette object at ...>
```

blow_out(*location=None*)

Force any remaining liquid to dispense, by moving this pipette's plunger to the calibrated *blow_out* position



Parameters: location (**Placeable** or tuple(**Placeable**, **Vector**)) – The **Placeable** (**Well**) to perform the blow_out. Can also be a tuple with first item **Placeable**, second item relative **Vector**

Returns:

Return This instance of [Pipette](#).
type:

Examples

```
>>> p200 = instruments.Pipette(name='p200', axis='a', max_volume=200)
>>> p200.aspirate(50).dispense().blow_out()
<opentrons.instruments.pipette.Pipette object at ...>
```

calibrate(*position*)

Calibrate a saved plunger position to the robot's current position

Notes

This will only work if the API is connected to a robot

Parameters: position (*str*) – Either "top", "bottom", "blow_out", or "drop_tip"

Returns:

Return This instance of [Pipette](#).
type:

Examples

```
>>> robot = Robot()
>>> p200 = instruments.Pipette(axis='a')
>>> robot.move_plunger(**{'a': 10})
>>> # save plunger 'top' to coordinate 10
>>> p200.calibrate('top')
<opentrons.instruments.pipette.Pipette object at ...>
```

calibrate_position(*location*, *current=None*)

Save the position of a **Placeable** (usually a **Container**) relative to this pipette.

Notes

The saved position will be persisted under this pipette's *name* and *axis* (see [Pipette](#))

Parameters: location (tuple(**Placeable**, **Vector**)) – A tuple with first item **Placeable**, second item relative **Vector**
current (**Vector**) – The coordinate to save this container to (Default: robot current position)

Returns:

Return This instance of [Pipette](#).
type:

Examples

```
>>> robot.reset()
<opentrons.robot.robot.Robot object at ...>
>>> tiprack = containers.load('tiprack-200ul', 'A1')
```



delay(seconds=0, minutes=0)

Parameters: seconds (*float*) – The number of seconds to freeeze in place.

dispense(volume=None, location=None, rate=1.0)

Dispense a volume of liquid (in microliters/uL) using this pipette

Notes

If no *location* is passed, the pipette will dispense from it's current position. If no *volume* is passed, *dispense* will default to it's *current_volume*

Parameters: volume (*int or float*) – The number of microliters to dispense (Default: self.current_volume)
location (**Placeable** or tuple(**Placeable**, **Vector**)) – The **Placeable** (**Well**) to perform the dispense. Can also be a tuple with first item **Placeable**, second item relative **Vector**
rate (*float*) – Set plunger speed for this dispense, where speed = rate * dispense_speed (see [set_speed\(\)](#))

Returns:

Return type: This instance of **Pipette**.

Examples

```
>>> p200 = instruments.Pipette(name='p200', axis='a', max_volume=200)
>>> # fill the pipette with liquid (200uL)
>>> p200.aspirate(plate[0])
<opentrons.instruments.pipette.Pipette object at ...>

>>> # dispense 50uL to a Well
>>> p200.dispense(50, plate[0])
<opentrons.instruments.pipette.Pipette object at ...>

>>> # dispense 50uL to the center of a well
>>> relative_vector = plate[1].center()
>>> p200.dispense(50, (plate[1], relative_vector))
<opentrons.instruments.pipette.Pipette object at ...>

>>> # dispense 20uL in place, at half the speed
>>> p200.dispense(20, rate=0.5)
<opentrons.instruments.pipette.Pipette object at ...>

>>> # dispense the pipette's remaining volume (80uL) to a Well
>>> p200.dispense(plate[2])
<opentrons.instruments.pipette.Pipette object at ...>
```

drop_tip(location=None, home_after=True)

Drop the pipette's current tip

Notes

If no location is passed, the pipette defaults to its *trash_container* (see [Pipette](#))



Returns:

Return type: This instance of [Pipette](#).

Examples

```
>>> robot.reset()
<opentrons.robot.robot.Robot object at ...>
>>> tiprack = containers.load('tiprack-200ul', 'A1')
>>> trash = containers.load('point', 'A1')
>>> p200 = instruments.Pipette(axis='a', trash_container=trash)
>>> p200.pick_up_tip(tiprack[0])
<opentrons.instruments.pipette.Pipette object at ...>
>>> # drops the tip in the trash
>>> p200.drop_tip()
<opentrons.instruments.pipette.Pipette object at ...>
>>> p200.pick_up_tip(tiprack[1])
<opentrons.instruments.pipette.Pipette object at ...>
>>> # drops the tip back at its tip rack
>>> p200.drop_tip(tiprack[1])
<opentrons.instruments.pipette.Pipette object at ...>
```

home()

Home the pipette's plunger axis during a protocol run

Notes

Pipette.home() homes the *Robot*

Returns:

Return type: This instance of [Pipette](#).

Examples

```
>>> p200 = instruments.Pipette(axis='a')
>>> p200.home()
<opentrons.instruments.pipette.Pipette object at ...>
```

mix(repetitions=1, volume=None, location=None, rate=1.0)

Mix a volume of liquid (in microliters/uL) using this pipette

Notes

If no *location* is passed, the pipette will mix from it's current position. If no *volume* is passed, *mix* will default to it's *max_volume*

Parameters:

- repetitions* (*int*) – How many times the pipette should mix (Default: 1)
- volume* (*int or float*) – The number of microliters to mix (Default: self.max_volume)
- location* ([Placeable](#) or tuple([Placeable](#), [Vector](#))) – The [Placeable](#) ([Well](#)) to perform the mix. Can also be a tuple with first item [Placeable](#), second item relative [Vector](#)
- rate* (*float*) – Set plunger speed for this mix, where speed = rate * (aspirate_speed or dispense_speed) (see [set_speed\(\)](#))

Returns:



examples

```
>>> p200 = instruments.Pipette(name='p200', axis='a', max_volume=200)

>>> # mix 50uL in a Well, three times
>>> p200.mix(3, 50, plate[0])
<opentrons.instruments.pipette.Pipette object at ...>

>>> # mix 3x with the pipette's max volume, from current position
>>> p200.mix(3)
<opentrons.instruments.pipette.Pipette object at ...>
```

`move_to(location, strategy='arc')`

Move this [Pipette](#) to a [Placeable](#) on the [Deck](#)

Notes

Until obstacle-avoidance algorithms are in place, [Robot](#) and [Pipette move_to\(\)](#) use either an “arc” or “direct”

Parameters: location ([Placeable](#) or tuple([Placeable](#), [Vector](#))) – The destination to arrive at
 strategy (“arc” or “direct”) – “arc” strategies (default) will pick the head up on Z axis, then over to the XY destination, then finally down to the Z destination. “direct” strategies will simply move in a straight line from the current position

Returns:

Return This instance of [Pipette](#).
 type:

`pick_up_tip(location=None, presses=3)`

Pick up a tip for the Pipette to run liquid-handling commands with

Notes

A tip can be manually set by passing a *location*. If no location is passed, the Pipette will pick up the next available tip in it's *tip_racks* list (see [Pipette](#))

Parameters: location ([Placeable](#) or tuple([Placeable](#), [Vector](#))) – The [Placeable \(Well\)](#) to perform the pick_up_tip. Can also be a tuple with first item [Placeable](#), second item relative [Vector](#)

Returns:

Return This instance of [Pipette](#).
 type:

Examples

```
>>> robot.reset()
<opentrons.robot.robot.Robot object at ...>
>>> tiprack = containers.load('tiprack-200ul', 'A1')
>>> p200 = instruments.Pipette(axis='a', tip_racks=[tiprack])
>>> p200.pick_up_tip(tiprack[0])
<opentrons.instruments.pipette.Pipette object at ...>
>>> p200.return_tip()
<opentrons.instruments.pipette.Pipette object at ...>
>>> # `pick_up_tip` will automatically go to tiprack[1]
>>> p200.pick_up_tip()
```



```
return_tip(home_alter=True)
```

Drop the pipette's current tip to it's originating tip rack

Notes

This method requires one or more tip-rack [Container](#) to be in this Pipette's *tip_racks* list (see [Pipette](#))

Returns:

Return type: This instance of [Pipette](#).

Examples

```
>>> tiprack = containers.load('tiprack-200ul', 'A1')
>>> p200 = instruments.Pipette(axis='a',
...     tip_racks=[tiprack], max_volume=200, name='p200')
>>> p200.pick_up_tip()
<opentrons.instruments.pipette.Pipette object at ...>
>>> p200.aspirate(50, plate[0])
<opentrons.instruments.pipette.Pipette object at ...>
>>> p200.dispense(plate[1])
<opentrons.instruments.pipette.Pipette object at ...>
>>> p200.return_tip()
<opentrons.instruments.pipette.Pipette object at ...>
```

set_speed(**kwargs)

Set the speed (mm/minute) the [Pipette](#) plunger will move during [aspirate\(\)](#) and [dispense\(\)](#)

Parameters: *kwargs (Dict)* – A dictionary who's keys are either "aspirate" or "dispense", and who's values are int or float
(Example: `{"aspirate": 300}`)

touch_tip(location=None, radius=1.0)

Touch the [Pipette](#) tip to the sides of a well, with the intent of removing left-over droplets

Notes

If no *location* is passed, the pipette will touch_tip from it's current position.

Parameters: *location (Placeable or tuple(Placeable, Vector))* – The [Placeable \(Well\)](#) to perform the touch_tip. Can also be a tuple with first item [Placeable](#), second item relative [Vector](#)
radius (float) – Radius is a floating point number between 0.0 and 1.0, describing the percentage of a well's radius. When radius=1.0, [touch_tip\(\)](#) will move to 100% of the wells radius. When radius=0.5, [touch_tip\(\)](#) will move to 50% of the wells radius.

Returns:

Return type: This instance of [Pipette](#).

Examples

```
>>> p200 = instruments.Pipette(name='p200', axis='a', max_volume=200)
>>> p200.aspirate(50, plate[0])
<opentrons.instruments.pipette.Pipette object at ...>
>>> p200.dispense(plate[1]).touch_tip()
<opentrons.instruments.pipette.Pipette object at ...>
```



Sign Up For Our Newsletter

Email**

SUBMIT

© OPENTRONS 2025